
BOLDQC


info@neuroinfo.org

Jun 09, 2023

CONTENTS:

| | | |
|----------|---|-----------|
| 1 | XNAT user documentation | 2 |
| 1.1 | Tagging your scans | 2 |
| 1.2 | Running the pipeline | 3 |
| 1.3 | Understanding the report page | 5 |
| 2 | Developer documentation | 16 |
| 2.1 | Installation | 16 |
| 2.2 | XNAT Installation | 17 |
| 3 | Indices and tables | 18 |

BOLDQC is a fMRI quality control pipeline. Working closely with neuroimaging experts, we designed a simple ergonomic user interface for the [XNAT](#) informatics and data management platform that allows users to quickly assess image quality and use those insights to address issues within the data acquisition workflow.



CHAPTER ONE

XNAT USER DOCUMENTATION

1.1 Tagging your scans

For BOLDQC to discover **BOLD** scans to process, you need to set the scan *type* in **XNAT**

| Type | Example series |
|------|----------------------|
| BOLD | ABCD_fMRI_rest_noPMU |

The image below displays an MR Session report page with populated types

PROJECT: TestProject01 > SUBJECT: 191118_SSBC_pilot > 191118_SSBC_pilot

MR Session: 191118_SSBC_pilot

Details Projects

Accession #: XNAT_E00004
Date Added: 2021-06-22 14:28:57 (admin)
Date: 2019-11-18
Time: 15:34:09
Operator: Mair
Scanner Name: AWP67056
Scanner Type: SIEMENS Prisma fit
Acquisition Site: CBS Neuroimaging

Subject: 191118_SSBC_pilot

Scans

| Bulk Actions: | | Type | Files | Note |
|--------------------------|------|-----------------------------|-----------------------|----------|
| <input type="checkbox"/> | Scan | ABCD_T1w_MPR_vNav_setter | 27.9 MB in 144 files | MOVE_001 |
| <input type="checkbox"/> | 12 | ABCD_T1w_MPR_vNav_setter | 42.9 MB in 176 files | ANAT_001 |
| <input type="checkbox"/> | 14 | ABCD_T1w_MPR_vNav | 441.6 MB in 383 files | |
| <input type="checkbox"/> | 17 | BOLD | | |
| <input type="checkbox"/> | 19 | ABCD_dMRI_Distortion Map_PA | 18.2 MB in 81 files | |

Total: 530.6 MB in 784 files

Actions

- Edit
- View
- Download
- Email
- Manage Files
- Delete
- Run Containers

1.2 Running the pipeline

To run the BOLDQC pipeline, use the Run Containers > boldqc-session button located within the Actions box on the MR Session report page

Note: If you don't see the Run Containers menu, please refer to [Setting up the container](#).

The screenshot shows the BOLDQC interface for an MR session. The top navigation bar includes 'Browse', 'New', 'Upload', 'Administer', 'Tools', and 'Help'. The user is logged in as 'admin'. The main content area shows the 'PROJECT: TestProject01 > SUBJECT: 191118_SSBC_pilot > 191118_SSBC_pilot' path. Below this, the 'MR Session: 191118_SSBC_pilot' is displayed. The 'Details' tab is selected, showing subject information: Accession #: XNAT_E00004, Subject: 191118_SSBC_pilot, Date Added: 2021-06-22 14:28:57 (admin), Gender: -, Date: 2019-11-18, Handedness: -, Time: 15:34:09, Age: --, Operator: Mair, Scanner Name: AWP67056, Scanner Type: SIEMENS Prisma_fit, and Acquisition Site: CBS Neuroimaging. To the right, an 'Actions' menu is open, with 'boldqc-session' highlighted under 'Run Containers'. The 'Scans' section below lists four scans with their details:

| | Scan | Type | Series Desc | Usability | Files | Note |
|--------------------------|------|-----------------------------|-----------------------------|-----------|-----------------------|----------|
| <input type="checkbox"/> | 12 | ABCD_T1w_MPR_vNav.setter | ABCD_T1w_MPR_vNav.setter | usable | 27.9 MB in 144 files | MOVE_001 |
| <input type="checkbox"/> | 14 | ABCD_T1w_MPR_vNav | ABCD_T1w_MPR_vNav | usable | 42.9 MB in 176 files | ANAT_001 |
| <input type="checkbox"/> | 17 | BOLD | ABCD_fMRI_rest | usable | 441.6 MB in 383 files | |
| <input type="checkbox"/> | 19 | ABCD_dMRI_Distortion Map_PA | ABCD_dMRI_Distortion Map_PA | usable | 18.2 MB in 81 files | |

Total: 530.6 MB in 784 files

This should bring up a small form with any configurable settings. Continue reading for a description of each setting

Set Container Launch Values

Please specify settings for this container.

session 191118_SSBC_pilot
session-label 191118_SSBC_pilot
run 1
Run number (Required)

Cancel Run Container

Manage Files
Delete
Run Containers

Scans

| Bulk Actions: | Download | Scan | Type | Series Desc | Usability | Files | Note |
|--------------------------|----------|------|----------------------|----------------------|-----------|-----------------------|----------|
| <input type="checkbox"/> | | 12 | ABCD_T1w_MPR_vNav | ABCD_T1w_MPR_vNav | usable | 27.9 MB in 144 files | MOVE_001 |
| <input type="checkbox"/> | | 14 | ABCD_T1w_MPR_vNav | ABCD_T1w_MPR_vNav | usable | 42.9 MB in 176 files | ANAT_001 |
| <input type="checkbox"/> | | 17 | BOLD | ABCD_fMRI_rest | usable | 441.6 MB in 383 files | |
| <input type="checkbox"/> | | 19 | ABCD_dMRI_Distortion | ABCD_dMRI_Distortion | usable | 18.2 MB in 81 files | Map_PA |

Total: 530.6 MB in 784 files

1.2.1 run

Note: The run number is *not* the scan number. If scan 17 was the first BOLD scan within the MR session, the run number would be 1, the second BOLD scan would be run 2, and so on.

This should be set to the integer value of the BOLD scan you want to process.

1.3 Understanding the report page

The following section will break down each section of the BOLDQC report page.

The screenshot shows the BOLDQC report page within an XNAT browser interface. The top navigation bar includes 'Browse', 'New', 'Upload', 'Administer', 'Tools', 'Help', 'Advanced', 'Search', and 'Logout'. The main content area displays the following sections:


- Summary**: Shows MR Session details (191118_SSBC_pilot), processing date (2021-10-28), and BOLD scan count (17).
- Parameters**: Lists Num Volumes (383), Num Voxels (91714), Mask Threshold (150), and Skip (4).
- QC Metrics**: Provides numerical values for Mean (247.036), StdDev (6.807), Slice SNR (229.1), Voxel SNR (42.571), Mean Rel Motion (0.104), Max Rel Motion (0.219), Mean Abs Motion (0.317), Max Abs Motion (0.467), Movements (>.1mm) (210), and Movements (>.5mm) (0).
- Files**: A list of generated files including SNR Image, Mean Image, StdDev Image, Mask Image, Mean Slice Image, Motion Data, Slope Image, Auto QC Report, and Slice Report.
- Images**: A tabbed section containing:
 - Signal-to-Noise Ratio (SNR)**: A grid of brain slices showing SNR across the volume.
 - Mean**: A grid of brain slices showing mean intensity across the volume.
 - Standard Deviation (StdDev)**: A grid of brain slices showing standard deviation across the volume.
 - Slope**: A grid of brain slices showing slope across the volume.
 - Mean Slice Intensity Over Time**: A line graph showing mean slice intensity over 350 volumes.
 - Motion Over Time**: Two line graphs showing rotations (pitch, roll, yaw) and translations (x, y, z) over 350 volumes.
 - Mask**: A grid of brain slices showing the segmentation mask.

1.3.1 Left pane

The left pane is broken up into several distinct sections. Each section will be described below.

Summary

The **Summary** pane orients the user to what MR Session they're currently looking at and various processing details



| Key | Description |
|------------|---------------------|
| MR Session | MR Session label |
| Processed | Processing date |
| BOLD scan | Processed BOLD scan |

Parameters

The Parameters pane displays fine-grained scan information

The screenshot shows the XNAT BOLDQC interface for a BOLD scan. The top navigation bar includes 'Browse', 'New', 'Upload', 'Administer', 'Tools', and 'Help'. The user is logged in as 'admin'.

The main content area displays the following information:


- Parameters:**
 - Num Volumes: 383
 - Num Voxels: 91714
 - Mask Threshold: 150
 - Skip: 4
- QC Metrics:**

| | |
|--------------------|---------|
| Mean: | 247.036 |
| StDev: | 6.807 |
| Slice SNR: | 229.1 |
| Voxel SNR: | 42.571 |
| Mean Rel Motion: | 0.104 |
| Max Rel Motion: | 0.219 |
| Mean Abs Motion: | 0.317 |
| Max Abs Motion: | 0.467 |
| Movements (>1mm): | 210 |
| Movements (>.5mm): | 0 |
- Files:**
 - SNR Image
 - Mean Image
 - StDev Image
 - Mask Image
 - Mean Slice Image
 - Motion Data
 - Slope Image
 - Auto QC Report
 - Slice Report
- Visualizations:**
 - Mean: A grid of 2D brain slices showing mean intensity across time points.
 - Standard Deviation (StDev): A grid of 2D brain slices showing standard deviation across time points.
 - Slope: A grid of 2D brain slices showing slope across time points.
 - Mean Slice Intensity Over Time: A line graph showing mean slice intensity over 383 volumes.
 - Motion Over Time: Two line graphs showing motion over time for rotations (pitch, roll, yaw) and translations (x, y, z) in mm.
 - Mask: A grid of 2D brain slices showing the mask used for analysis.

| Key | Description |
|----------------|---|
| Num Volumes | Number of time points |
| Num Voxels | Number of voxels included in the analysis |
| Mask Threshold | Masking threshold |
| Skip | Number of initial time points discarded |

QC Metrics

The QC Metrics pane displays quality control metrics computed *over the entire volume*



| Metric | Description |
|-------------------|--|
| Mean | Mean signal intensity |
| StDev | Mean voxel Standard deviation |
| Slice SNR | Mean slice-based SNR (sensitive to motion) |
| Voxel SNR | Mean voxel SNR |
| Mean Rel Motion | Mean relative translations in 3D (mm) |
| Max Rel Motion | Maximum relative motion (mm) |
| Mean Abs Motion | Mean absolute motion in 3D (mm) |
| Max Abs Motion | Maximum absolute motion in 3D (mm) |
| Movements (>.1mm) | Number of relative translations in 3D > .1mm |
| Movements (>.5mm) | Number of relative translations in 3D > .5mm |

Files

The **Files** pane contains the most commonly requested files. Clicking on any of these files will display that file in the browser


| File | Description |
|------------------|---|
| SNR Image | BOLD signal SNR image, axial |
| Mean Image | BOLD signal mean image, axial |
| StDev Image | BOLD signal standard deviation image, axial |
| Mask Image | Masked image |
| Mean Slice Image | BOLD signal mean slice intensity plot |
| Motion Data | Motion (translations and rotations) plot |
| Slope Image | BOLD signal slope image, axial |
| Auto QC Report | Automated QC report |
| Slice Report | Individual slice QC report |

1.3.2 Tabs


To the right of the *left pane* left pane you'll find a tab container.

Images

The Images tab displays a zoomed out view of the SNR, Mean, Standard Deviation, Slope, MEan Slice Intensity, Motion, and Mask images



Clicking on any of these images will display a larger version of the image



Automated QC

The **Automated QC** tab displays a complete list of BOLDQC metrics.

PROJECT: TestProject01 > SUBJECT: 191118_SSBC_pilot > SESSION: 191118_SSBC_pilot > 191118_SSBC_pilot_BOLD_17_EQC

| Action | > | |
|--------------------|-------------------|--|
| Summary | | |
| MR Session | 191118_SSBC_pilot | |
| Processed | 2021-10-28 | |
| BOLD scan | 17 | |
| Parameters | | |
| Num Volumes: | 383 | |
| Num Voxels: | 91714 | |
| Mask Threshold: | 150 | |
| Skip: | 4 | |
| QC Metrics | | |
| Mean: | 247.036 | |
| StdDev: | 6.807 | |
| Slice SNR: | 229.1 | |
| Voxel SNR: | 42.571 | |
| Mean Rel Motion: | 0.104 | |
| Max Rel Motion: | 0.219 | |
| Mean Abs Motion: | 0.317 | |
| Max Abs Motion: | 0.467 | |
| Movements (>.1mm): | 210 | |
| Movements (>.5mm): | 0 | |
| Files | | |
| SNR Image | | |
| Mean Image | | |
| StdDev Image | | |
| Mask Image | | |
| Mean Slice Image | | |
| Motion Data | | |
| Slope Image | | |
| Auto QC Report | | |
| Slice Report | | |

| Images | Automated QC | Manual QC | All Stored Files |
|--------|----------------|------------|------------------|
| # | Metric | Value | |
| 1 | N_Vols | 383 | |
| 2 | Skip | 4 | |
| 3 | qc_N_Tps | 379 | |
| 4 | qc_thresh | 150.000000 | |
| 5 | qc_nVox | 91714 | |
| 6 | qc_Mean | 247.035605 | |
| 7 | qc_Max | 315.690000 | |
| 8 | qc_Min | 94.540000 | |
| 9 | qc_Stddev | 6.806792 | |
| 10 | qc_sSNR | 229.100000 | |
| 11 | qc_vSNR | 42.570848 | |
| 12 | qc_slope | 0.000476 | |
| 13 | mot_N_Tps | 379 | |
| 14 | mot_rel_x_mean | 0.005368 | |
| 15 | mot_rel_x_sd | 0.004458 | |
| 16 | mot_rel_x_max | 0.026579 | |
| 17 | mot_rel_x_1mm | 0 | |
| 18 | mot_rel_x_5mm | 0 | |
| 19 | mot_rel_y_mean | 0.092484 | |
| 20 | mot_rel_y_sd | 0.049973 | |
| 21 | mot_rel_y_max | 0.212603 | |
| 22 | mot_rel_y_1mm | 176 | |
| 23 | mot_rel_y_5mm | 0 | |
| 24 | mot_rel_z_mean | 0.036497 | |
| 25 | mot_rel_z_sd | 0.024515 | |

Manual QC

The Manual QC tab contains a form allowing a quality control technician to record additional observations, comments, and assign a final PASS, WARN, or FAIL grade to the scan

The screenshot shows the XNAT web interface for the BOLDQC application. The URL is <https://ncftest01.rc.fas.harvard.edu>. The top navigation bar includes links for Browse, New, Upload, Administer, Tools, Help, Advanced, and Logout. The user is logged in as 'admin'.

The main content area shows a detailed QC report for a session named '191118_SSBC_pilot'. The report includes sections for Summary, Parameters, QC Metrics, and Files. The 'Manual QC' tab is selected, showing sections for Motion, Head Coverage, Unexpected Inhomogeneity, and Ghosting/Blurring - Affects Brain. Each section contains a list of options for the user to select.

Summary

- MR Session: 191118_SSBC_pilot
- Processed: 2021-10-28
- BOLD scan: 17

Parameters

| | |
|-----------------|-------|
| Num Volumes: | 383 |
| Num Voxels: | 91714 |
| Mask Threshold: | 150 |
| Skip: | 4 |

QC Metrics

| | |
|--------------------|---------|
| Mean: | 247.036 |
| StdDev: | 6.807 |
| Slice SNR: | 229.1 |
| Voxel SNR: | 42.571 |
| Mean Rel Motion: | 0.104 |
| Max Rel Motion: | 0.219 |
| Mean Abs Motion: | 0.317 |
| Max Abs Motion: | 0.467 |
| Movements (>.1mm): | 210 |
| Movements (>.5mm): | 0 |

Files

- SNR Image
- Mean Image
- StdDev Image
- Mask Image
- Mean Slice Image
- Motion Data
- Slope Image
- Auto QC Report
- Slice Report

Motion

- Good
- Questionable
- Bad

Head Coverage

- Good - full brain coverage with ≥ 1 slice buffer
- Questionable - slight clipping or necessary clipping
- Bad - brain clipped unnecessarily or severely

Unexpected Inhomogeneity

- None
- Expected
- Unexpected - moderate
- Unexpected - severe

Ghosting/Blurring - Affects Brain

- None
- Mild
- Severe

POWERED BY

All Stored Files

The All Stored Files tab contains a list of *every file* stored by BOLDQC

The screenshot shows a web browser window for the XNAT platform. The URL is <https://ncf01.rc.fas.harvard.edu>. The top navigation bar includes links for Browse, New, Upload, Administer, Tools, Help, Advanced, and Search. A user is logged in as 'admin'.

The main content area shows a report for a BOLD scan session. The session details are:

- MR Session:** 191118_SSBC_pilot
- Processed:** 2021-10-28
- BOLD scan:** 17

Parameters:

| | |
|-----------------|-------|
| Num Volumes: | 383 |
| Num Voxels: | 91714 |
| Mask Threshold: | 150 |
| Skip: | 4 |

QC Metrics:

| | |
|--------------------|---------|
| Mean: | 247.036 |
| StdDev: | 6.807 |
| Slice SNR: | 229.1 |
| Voxel SNR: | 42.571 |
| Mean Rel Motion: | 0.104 |
| Max Rel Motion: | 0.219 |
| Mean Abs Motion: | 0.317 |
| Max Abs Motion: | 0.467 |
| Movements (>.1mm): | 210 |
| Movements (>.5mm): | 0 |

Files:

| |
|------------------|
| SNR Image |
| Mean Image |
| StdDev Image |
| Mask Image |
| Mean Slice Image |
| Motion Data |
| Slope Image |
| Auto QC Report |
| Slice Report |

All Stored Files:

| # | File | Description | Type |
|----|---|----------------------------|--------------|
| 1 | 191118_SSBC_pilot_BOLD_17_EQC_auto_report.txt | Automated QC Report | text/plain |
| 2 | 191118_SSBC_pilot_BOLD_17_EQC_mask_thumbnail.png | Mask Image | image/png |
| 3 | 191118_SSBC_pilot_BOLD_17_EQC_mean.nii.gz | Mask NIFTI | image/nifti1 |
| 4 | 191118_SSBC_pilot_BOLD_17_EQC_mean_thumbnail.png | Mean Image | image/png |
| 5 | 191118_SSBC_pilot_BOLD_17_EQC_mean.nii.gz | Mean NIFTI | image/nifti1 |
| 6 | 191118_SSBC_pilot_BOLD_17_EQC_mean_slice.txt | Mean Slice Intensity Data | text/plain |
| 7 | 191118_SSBC_pilot_BOLD_17_EQC_mean_slice.png | Mean Slice Intensity Image | text/html |
| 8 | 191118_SSBC_pilot_BOLD_17_EQC_motion.png | Motion Image | image/png |
| 9 | 191118_SSBC_pilot_BOLD_17_EQC_slice_report.txt | Slice Report | text/plain |
| 10 | 191118_SSBC_pilot_BOLD_17_EQC_slope_thumbnail.png | Slope Image | text/png |
| 11 | 191118_SSBC_pilot_BOLD_17_EQC_slope.nii.gz | Slope NIFTI | image/nifti1 |
| 12 | 191118_SSBC_pilot_BOLD_17_EQC_snrr_thumbnail.png | SNR Image | image/png |
| 13 | 191118_SSBC_pilot_BOLD_17_EQC_snrr.nii.gz | SNR NIFTI | image/nifti1 |
| 14 | 191118_SSBC_pilot_BOLD_17_EQC_stdev_thumbnail.png | StdDev Image | image/png |
| 15 | 191118_SSBC_pilot_BOLD_17_EQC_stdev.nii.gz | StdDev NIFTI | image/nifti1 |

POWERED BY

Note: Clicking on a file within the All Stored Files tab will download the file.

| File | Description |
|---------------------------|-----------------------------------|
| *_EQC_auto_report.txt | Automated QC report |
| *_EQC_mask_thumbnail.png | Mask snapshot image |
| *_EQC_mask.nii.gz | Mask NIFTI |
| *_EQC_mean_thumbnail.png | Mean snapshot image |
| *_EQC_mean.nii.gz | Mean NIFTI |
| *_EQC_mean_slice.txt | Mean slice intensity data |
| *_EQC_mean_slice.png | Mean slice intensity plot |
| *_EQC_motion.png | Motion plot |
| *_EQC_slice_report.txt | Slice report |
| *_EQC_slope_thumbnail.png | Slope snapshot image |
| *_EQC_slope.nii.gz | Slope NIFTI |
| *_EQC_snr_thumbnail.png | SNR snapshot image |
| *_EQC_snr.nii.gz | SNR NIFTI |
| *_EQC_stdev_thumbnail.png | Standard deviation snapshot image |
| *_EQC_stdev.nii.gz | Standard deviation image |

DEVELOPER DOCUMENTATION

2.1 Installation

At the moment, the only supported way to install BOLDQC is *within a container*.

2.1.1 downloading a container

There are prebuilt versions of BOLDQC on Docker Hub. You can pull the latest version by running

```
docker pull neuroinformatics/boldqc
```

or you can pull a specific version e.g., `0.1.0` by running

```
docker pull neuroinformatics/boldqc:0.1.0
```

2.1.2 building a container

To build BOLDQC as a container, grab the latest `Dockerfile` from the repository and run

```
docker build -t boldqc:latest - < Dockerfile
```

Now you can run `boldQC.py`—which is the default `ENTRYPOINT`—using `docker run`

```
docker run boldqc:latest --help
```

Note: You can also convert the BOLDQC Docker image into a Singularity image, however to run `boldQC.py` you'll need to supply `--pwd /sw/apps/boldqc`

```
singularity run --pwd /sw/apps/boldqc boldqc.sif --help
```

2.2 XNAT Installation

The following section will describe how to build and configure BOLDQC as a [XNAT](#) plugin.

2.2.1 building the plugin

Clone the `xnat-1.8` branch from the github.com/harvard-nrg/boldqc repository

```
git clone -b xnat-1.8 --single-branch https://github.com/harvard-nrg/boldqc
```

Change into the repository directory and compile the plugin using [Gradle](#)

```
./gradlew jar
```

Once the plugin has been compiled, move the resulting `.jar` into your XNAT plugins directory

```
mv ./build/libs/boldqc-plugin-1.0.0.jar ${XNAT_HOME}/plugins/
```

**CHAPTER
THREE**

INDICES AND TABLES

- genindex
- modindex
- search